

CeeMax™ 901ss Table Top Slit Lamp

Slit Illumination

Type	Galilean-Type
Magnification Ratio	10X, 16X, 25X
Objective	F: 108mm
Eyepieces	12.5X
Angle between eyepieces	13° Converging optics
Pupillary Adjustment	55mm~78.5mm
Diopter	-6 to+6 Diopter
Slit Width (mm)	Continuously variable from 0 to 14 (at 14mm, slit becomes a circle)
Slit Length (mm)	Continuously variable from 1 to 14
Aperture Diameters (mm)	14, 8, 5, 3, 0.5
Slit Angle	0°~180° Continuously variable
Filters	Heat absorbing, UV, Red-free, Cobalt Blue
Lamp	10W LED
Illumination Type	From Bottom (Zeiss type)

Chin Rest

Vertical movement	70mm
Fixation Target	Green LED

Base

Longitudinal movement	90mm
Lateral movement	100mm
Fine Base movement	15mm
Vertical movement	30m

Demension & Weight

Package Dimension 730mmx430mmx470mm
Total Weight 24Kg (53 lb)

Power Supply

Input 110V/230V AC
Frequency 60Hz/50Hz



- Ready for digital photography and video
- Ready for Goldmann tonometer
- One year factory warranty (contact CeeMax™ Ophthalmic)

CeeMax™ Ophthalmic

17100 Gillette Ave. Unit 113, Irvine, CA 92614 USA

Tel: 949.885.8875

Email: info@ceemaxusa.com