

# CeeMax™ 902ss Table Top Slit Lamp

## Slit Illumination

Type	Galilean-Type
Magnification Ratio	6X, 10X, 16X, 25X, 40X
Objective	F: 108mm
Eyepieces	12.5X
Angle between eyepieces	13° Converging optics
Pupillary Adjustment	55mm~78.5mm
Diopter	-6 to+6 Diopter
Slit Width (mm)	Continuously variable from 0 to 14 (at 14mm, slit becomes a circle)
Slit Length (mm)	Continuously variable from 1 to 14
Aperture Diameters (mm)	14, 8, 5, 3, 0.5
Slit Angle	0°~180° Continuously variable
Filters	Heat absorbing, UV, Red-free, Cobalt Blue
Lamp	10W LED
Illumination Type	From Bottom (Zeiss type)

## Chin Rest

Vertical movement	70mm
Fixation Target	Green LED

## Base

Longitudinal movement	90mm
Lateral movement	100mm
Fine Base movement	15mm
Vertical movement	30m

## Demension & Weight

Package Dimension 730mmx430mmx470mm  
Total Weight 24Kg (53 lb)

## Power Supply

Input 110V/230V AC  
Frequency 60Hz/50Hz



- Ready for digital photography and video
- Ready for Goldmann tonometer
- One year factory warranty (contact CeeMax™ Ophthalmic)

**CeeMax™ Ophthalmic**

**17100 Gillette Ave. Unit 113, Irvine, CA 92614 USA**

**Tel: 949.885.8875**

**Email: [info@ceemaxusa.com](mailto:info@ceemaxusa.com)**